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FOREWORD

The material is intended for students in laboratory exercises in the subject Food Chemistry of the academic study programme Food Science and Nutrition and is a supplement to the existing study material "Abram and Zelenik-Blatnik, Exercises in food chemistry for students of food technology, 2002". The additional study material summarizes and at the same time builds on the theoretical starting points and extends the research approach of the laboratory exercises described in the aforementioned basic study material. The authors believe that a clear schematic presentation of the theoretical basis and practical implementations, as well as a comprehensive explanation of the processing of experimental data, which we have included in the additional study material, are a great help in understanding a relatively complex issue.

Students from an international environment (Erasmus exchange) participate in the Food Chemistry course at the same time as regular students. We have prepared an English version of the material to offer those students suitable material that is harmonized with the Slovenian one.

With this material, we would like to encourage students to effectively prepare and carry out laboratory exercises in the subject Food Chemistry, increase their safety when working in the laboratory, strengthen their understanding and independence in their work, and improve their key skills.

BASICS OF CHEMICAL CALCULATIONS

Concentration units

Mass fraction of the solute (w_{solute}) indicates the ratio between the mass of the solute (m_{solute}) and the total mass of the solution (m_{solution}). Unit: %.

$$w_{\text{solute}} = \frac{m_{\text{solute}}}{m_{\text{solution}}} \times 100 \%$$

For a solution prepared as a mixture of different liquids, the liquid solute content is expressed as a volume fraction (φ_{solute}). The volume fraction of the solute indicates the ratio between the volume of the pure liquid solute (V_{solute}) and the total volume of the solution (V_{solution}). Unit: %.

$$\varphi_{\text{solute}} = \frac{V_{\text{solute}}}{V_{\text{solution}}} \times 100 \%$$

Mass concentration (γ) indicates the mass of solute in a certain volume of solution. Unit: g/mL, mg/mL, g/L,

$$\gamma = \frac{m_{\text{solute}}}{V_{\text{solution}}}$$

Molar concentration (c) indicates the number of moles of solute in 1 L of solution. Unit: mol/L.

$$c = \frac{n_{\text{solute}}}{V_{\text{solution}}}$$

When we dilute a more concentrated solution (stock solution), we measure a certain volume (usually with a pipette) of the stock solution ($V_{p \text{ stock}}$) and add sufficient solvent to prepare the corresponding volume of the diluted solution ($V_{\text{diluted slution}}$). The solute concentration in a diluted solution is calculated:

$$\text{solute concentration in diluted solution} = \frac{\text{solute concentration in stock solution} \times V_{p \text{ stock}}}{V_{\text{diluted slution}}}$$

DETERMINATION OF GLUCOSE AND SUCROSE IN THE MIXTURE

Glucose

D(+)-glucose (grape sugar or dextrose) is an aldohexose with the molecular formula $C_6H_{12}O_6$.

Next to fructose, glucose is the most common monosaccharide in foods. It is present in a relatively high proportion in fruit, where its content increases with the degree of ripeness. As a monomeric unit, it is a building block of many disaccharides (e.g. sucrose, lactose, maltose), oligo- and polysaccharides (e.g. β -glucans, amylose and amylopectin (starch), cellulose), and compounds that are formed in the reaction between the hemiacetal $-OH$ group of glucose and a non-sugar molecule (called aglycone), e.g. with alcohol, a phenolic compound,

In the food industry, glucose is a raw material for the production of the polyol sorbitol (reduction of the carbonyl group), gluconic acid and lactone (oxidation of the carbonyl group), and for the production of ethanol (by fermentation).

In the presence of oxidizing agents such as the *Luff* reagent and the *Fehling* reagent (alkaline solution of Cu^{2+} ions), glucose (carbonyl group) is oxidized to carboxylic acid. Here, Cu^{2+} is reduced to Cu^+ . This is why we say that glucose is a reducing sugar.

In an alkaline environment (or enzyme-catalysed), the isomerization of ketose to aldose takes place. Under these conditions, even fructose (which, like ketose, is not a reducing sugar) acts as a "reducing sugar".

Sucrose

Sucrose (table sugar, cane sugar, beet sugar) is α -D-glucopyranosyl- β -D-fructofuranoside with the molecular formula $C_{12}H_{22}O_{11}$. D-glucopyranose and D-fructofuranose are linked by an α,β -(1 \rightarrow 2)-glycosidic bond. The reducing ends of both molecules are involved in a glycosidic bond; therefore, sucrose is not a reducing sugar (meaning it is less reactive).

Sucrose is the most abundant disaccharide in plants, where it plays a role in the transport of sugars and as an energy store. It accumulates in larger amounts in the vegetative parts of the plant, such as leaves and stems (sugar cane up to about 25%), and in fruit (apples, oranges, apricots, peaches up to about 8%). Some fruits (e.g. grapes, figs) contain small amounts of sucrose, as this is hydrolyzed during ripening. In vegetables, the sucrose content is considerable in sugar beet (up to about 20%) and onions (up to 10%). In potatoes, the sucrose content in the tubers is relatively high before full ripening and decreases during ripening.

The industrial source of sucrose is sugar beet and sugar cane. In some geographical areas, other sources are also used (e.g. dates, maple sap, ...).

Under acidic conditions or in the presence of the enzyme β -fructofuranosidase (invertase), the glycosidic bond in the sucrose molecule is cleaved (hydrolysed). When the glycosidic bond in the sucrose molecule is cleaved, a mixture of glucose and fructose is formed. The cleavage of the bonds in one mole of sucrose produces two moles of monosaccharides; in addition, the relative sweetness of glucose-fructose syrup is higher than that of sucrose. Therefore, the solution of the resulting mixture is sweeter than the sucrose solution. The name invert sugar for sucrose is due to the inversion: the specific rotation of linearly polarised light for sucrose is $+66.5^\circ$, and for the mixture of monosaccharides after the hydrolysis of sucrose it is -33.3° .

Sucrose is used as a sweetener in foods because it has a sweet taste (like most mono- and disaccharides). Due to the high number of hydrophilic $-OH$ groups in the molecule, sucrose (like other mono- and disaccharides and shorter oligosaccharides) is very soluble in water. This enables the preparation of very concentrated solutions (syrup, honey) with high osmolarity, which inhibits the growth of microorganisms. Sucrose in food can also act as a humectant.

In the diet, sucrose is a source of energy. The enzyme β -fructofuranosidase in the mucosa of the small intestine catalyses the cleavage of the glycosidic bond. The resulting monosaccharides are absorbed and transported into the blood.

Objective

Determine the mass of glucose (m_{glucose}) and the mass of sucrose (m_{sucrose}) in the mixture by the *Schoorl-Luff* method.


Principle of method

The determination of glucose is based on an oxidation-reduction reaction that takes place after adding *Luff* reagent (alkaline solution of copper citrate complex) when glucose (reducing sugar) reduces Cu^{2+} to Cu_2O in an alkaline medium.


The determination of sucrose (not a reducing sugar) in the mixture is based on the cleavage of the glycosidic bond (hydrolysis) and the determination of the resulting hydrolysis products (glucose and fructose) and the glucose already present in the sample by adding *Luff* reagent. Here, the isomerization of fructose to a reducing sugar takes place in an alkaline environment.

When a reducing sugar solution is added to *Luff* reagent, an oxidation-reduction reaction takes place, in which the Cu_2O is formed as a reddish-brown precipitate.

When iodine solution is added, the following oxidation-reduction reaction takes place in the excess amount of iodine:

We perform a back titration by determining the amount of remaining I_2 ($n_{\text{I}_2 \text{ remaining}}$; excess with respect to the amount of Cu_2O) by titration with $\text{Na}_2\text{S}_2\text{O}_3$ solution; this is the sample analysis.

The total amount of added I_2 ($n_{\text{I}_2 \text{ total}}$) is determined by the blank sample analysis, which is performed in the same way as the sample analysis, except that water is added to *Luff* reagent instead of reducing sugar.

The amount of reducing sugar or the amount of Cu_2O formed in the reaction between the reducing sugar and *Luff* reagent is determined as:


$$n_{\text{I}_2 \text{ reacted}} = n_{\text{I}_2 \text{ total}} - n_{\text{I}_2 \text{ remaining}}$$

Experimental procedure

- a) Diluted sample preparation
- b) Glucose determination
- c) Hydrolysis - after hydrolysis of the glycosidic bond in sucrose, the hydrolysate contains the hydrolysis products and the glucose already present in the sample
- d) Determination of all reducing sugars in hydrolysate (glucose already present in the sample + reducing sugars formed by hydrolysis)
- e) Blank sample analysis (water instead of the sugar solution)

a) Diluted sample preparation


- The sample in a 250 mL volumetric flask containing glucose and sucrose is filled to the mark with distilled water, the flask is closed and the contents are mixed - the volume of the diluted sample ($V_{\text{diluted sample}} = 250 \text{ mL}$).

b) Glucose determination

- Measure 25 mL of *Luff* reagent solution into a glass-stoppered Erlenmeyer flask and add 25 mL of the diluted sample.
- Place the flask on a hot plate and connect it to a water-cooled condenser; heat at boiling temperature for 10 min; remove from the condenser; cool to room temperature.
- After cooling, add: 50 mL of CH₃COOH solution (0.4 mol/L; measuring cylinder), 10 mL of iodine solution (0.1 mol/L; volumetric pipette) and 55 mL of HCl solution (0.75 mol/L; measuring cylinder); shake the contents (first slowly and carefully, then more intensively).
- Titrate with a standard aqueous solution of Na₂S₂O₃ (0.1 mol/L) to the color change from dark blue to light blue using 1 mL of starch solution (1%) as an indicator.
- At the endpoint of the titration, read the volume of the Na₂S₂O₃ solution, i.e. V_b.


SCHEME OF GLUCOSE DETERMINATION PROCEDURE

c) Hydrolysis

- Measure 50 mL of the diluted sample from the 250 mL flask into a 100 mL flask.
- Add 5 mL of HCl solution (36%; HCl is concentrated; care must be taken).
- Place the flask in a water bath at 70 °C for 10 min (acid hydrolysis of the glycosidic bond in the sucrose molecule is taking place); then cool the contents of the flask.
- Add phenolphthalein and neutralize with NaOH solution (15%).
- Fill the contents in the 100 mL flask (hydrolysate) to the mark with distilled water.
- The products of the cleavage of the glycosidic bond in sucrose (glucose and fructose) and glucose already present in the sample are in a 100 mL flask (in the hydrolysate).

SCHEME OF HYDROLYSIS PROCEDURE

d) Determination of all reducing sugars in hydrolysate

- Measure 25 mL of *Luff* reagent solution into a glass-stoppered Erlenmeyer flask and add 25 mL of the hydrolysate from a 100 mL flask.
- Proceed in the same way as for the determination of glucose (scheme p. 133).
- At the endpoint of titration, read the volume of the $\text{Na}_2\text{S}_2\text{O}_3$ solution, i.e. V_c .

e) Blank sample analysis

- Measure 25 mL of *Luff* reagent solution into a glass-stoppered Erlenmeyer flask and add 25 mL of distilled water.
- Proceed in the same way as for the determination of glucose (scheme p. 133).
- At the endpoint of titration, read the volume of the $\text{Na}_2\text{S}_2\text{O}_3$ solution, i.e. V_a .

Calculations

From titration data we calculate:

- the amount of $\text{Na}_2\text{S}_2\text{O}_3$ which is proportional to $n_{\text{I}_2 \text{ reacted}}$ in the determination of m_{glucose} in 25 mL of diluted sample ($n_{\text{Na}_2\text{S}_2\text{O}_3 \text{ glucose}}$):

$$n_{\text{Na}_2\text{S}_2\text{O}_3 \text{ glucose}} = (V_a - V_b) \times c_{\text{Na}_2\text{S}_2\text{O}_3}$$

- the amount of $\text{Na}_2\text{S}_2\text{O}_3$ which is proportional to $n_{\text{I}_2 \text{ reacted}}$ in the determination of $m_{\text{all reduc sugars}}$ in 25 mL of hydrolysate ($n_{\text{Na}_2\text{S}_2\text{O}_3 \text{ all reduc sugars}}$):

$$n_{\text{Na}_2\text{S}_2\text{O}_3 \text{ all reduc sugars}} = (V_a - V_c) \times c_{\text{Na}_2\text{S}_2\text{O}_3}$$

To calculate the mass of the sugar, we do not use the stoichiometric ratio between the amount of sugar and the amount of $n_{\text{Na}_2\text{S}_2\text{O}_3}$, but the experimentally determined ratio (Table 1 on p. 138).

Using the data in Table 1 on p. 138, we draw a graph and read:

- m_{glucose} in 25 mL of diluted sample at calculated $n_{\text{Na}_2\text{S}_2\text{O}_3 \text{ glucose}}$ and
- $m_{\text{all reduc sugars}}$ in 25 mL of hydrolysate at calculated $n_{\text{Na}_2\text{S}_2\text{O}_3 \text{ all reduc sugars}}$.

Table 1: Conversion of the $n_{\text{Na}_2\text{S}_2\text{O}_3}$ amount to the mass of reducing sugar in 25 mL of diluted sample or hydrolysate for the determination method described.

$n_{\text{Na}_2\text{S}_2\text{O}_3}$ (mmol)	$m_{\text{reduc sugar}}$ (mg)
0.1	2.4
0.2	4.8
0.3	7.2
0.4	9.7
0.5	12.2
0.6	14.7
0.7	17.2
0.8	19.8

Calculation of the mass of glucose in the sample

According to the experiment description: 25 mL of the diluted sample is measured from a 250 mL flask into a glass-stoppered Erlenmeyer flask, therefore:

$$m_{\text{glucose}} \text{ in 250 mL of diluted sample} = m_{\text{glucose}} \text{ in 25 mL of diluted sample} \times 10$$

$$m_{\text{glucose}} \text{ in the sample} = m_{\text{glucose}} \text{ in 250 mL of diluted sample}$$

Calculation of the mass of sucrose in the sample

According to the experiment description: 25 mL of hydrolysate is measured from a 100 mL flask containing hydrolysate, therefore:

$$\begin{aligned} m_{\text{all reduc sugars}} \text{ in 100 mL of hydrolysate} &= m_{\text{all reduc sugars}} \text{ in 25 mL of hydrolysate} \times 4 = \\ &= m_{\text{all reduc sugars}} \text{ in 50 mL of diluted sample measured into 100 mL flask} \end{aligned}$$

According to the experiment description: 50 mL of diluted sample is measured from a 250 mL flask, therefore:

$$m_{\text{all reduc sugars}} \text{ in 250 mL of diluted sample} = m_{\text{all reduc sugars}} \text{ in 50 mL of diluted sample} \times 5$$

$$m_{\text{all reduc sugars}} \text{ in the sample} = m_{\text{all reduc sugars}} \text{ in 250 mL of diluted sample}$$

We take into account the fact that water is released when two monosaccharides are joined to form a disaccharide, therefore: $M_{\text{sucrose}} = 2 \times 180 \text{ g/mol} - 18 \text{ g/mol}$

$$\text{mass of disaccharide} = \text{mass of two monosaccharides} \times 0.95$$

$$m_{\text{sucrose}} \text{ in the sample} = (m_{\text{all reduc sugars}} \text{ in the sample} - m_{\text{glucose}} \text{ in the sample}) \times 0.95$$

Report

Experimental data


$$V_a = \dots\dots\dots V_b = \dots\dots\dots V_c = \dots\dots\dots C_{\text{Na}_2\text{S}_2\text{O}_3} = \dots\dots\dots$$

Further data can be found in the description of the experimental procedure.

Calculations and graph

$$n_{\text{Na}_2\text{S}_2\text{O}_3 \text{ glucose}} =$$

$$n_{\text{Na}_2\text{S}_2\text{O}_3 \text{ all reduc sugars}} =$$

Dependence of reducing sugar mass on $\text{Na}_2\text{S}_2\text{O}_3$ amount

The mass of glucose in 25 mL of diluted sample read from the graph at calculated value of

$$n_{\text{Na}_2\text{S}_2\text{O}_3 \text{ glucose}} = \dots\dots\dots$$

$$m_{\text{glucose}} \text{ in 250 mL of diluted sample} =$$

$$m_{\text{glucose}} \text{ in the sample} =$$

The mass of all reducing sugars in 25 mL of hydrolysate read from the graph at calculated

$$\text{value of } n_{\text{Na}_2\text{S}_2\text{O}_3 \text{ all reduc sugars}} = \dots\dots\dots$$

$$m_{\text{all reduc sugars}} \text{ in 100 mL of hydrolysate} =$$

$$m_{\text{all reduc sugars}} \text{ in 250 mL of diluted sample} =$$

$$m_{\text{all reduc sugars}} \text{ in the sample} =$$

$$m_{\text{sucrose}} \text{ in the sample} =$$

Result

sample number

m_{glucose} in the sample =

m_{sucrose} in the sample =

DETERMINATION OF TRYPSIN INHIBITOR ACTIVITY IN SOY

Trypsin inhibitors

Enzyme inhibitors are various substances that reduce the activity of enzymes and occur naturally or are added to foods. From a nutritional point of view, inhibitors of proteolytic digestive enzymes, the proteases, are particularly important.

Protease inhibitors are found in plants (mainly in legumes, but also in cereals, potatoes, tomatoes) where they are supposed to fulfill various tasks - protecting the plant from endogenous proteases that are released when the plant tissue is damaged; protecting against proteases from microorganisms, insects, The highest content of these inhibitors was found in seeds, leaves and tubers. Protease inhibitors of microbial or animal origin are also known.

Protease inhibitors are often proteins or polypeptides. Depending on the catalytic type of protease they inhibit, protease inhibitors are classified into four groups.

The soybean serine protease inhibitors include:

- Kunitz-type inhibitors contain few cysteine residues and form two disulfide bridges. They are less thermally stable.
- Bowman-Birk type inhibitors contain several cysteine residues and form several disulfide bridges. They are more thermally stable. Compared with Kunitz-type inhibitors, they have lower molar mass.

Both types of serine protease inhibitors are potent trypsin inhibitors.

Trypsin is an endoprotease that catalyses the cleavage of the peptide bond in proteins and peptide substrates on the carboxyl side of the amino acids lysine or arginine. The inhibition of trypsin is the result of the formation of a very stable complex between the inhibitor and the enzyme.

In humans, the inhibition of trypsin leads to incomplete protein digestion, which reduces the bioavailability of amino acids from proteins. Inhibition of trypsin leads to excessive production and secretion of trypsin, which in extreme cases can lead to hypertrophy (enlargement) of the pancreas. In addition, the increased synthesis of pancreatic enzymes impairs synthesis of other proteins in the body, resulting in slower growth.

The activity of trypsin inhibitors is reduced by thermal treatment, as denaturation occurs at elevated temperatures due to disruption of the native structure. In the presence of reducing agents, where the disulfide bonds in the inhibitors are reduced, thermal treatment leads to even more effective denaturation.

Since the use of legumes, which are particularly rich in trypsin inhibitors (especially soy), is quite widespread in the diet, care must be taken to ensure that they are sufficiently thermally treated. The method, time and temperature of thermal treatment, the particle size and the moisture content of the seeds influence the extent to which the activity of the trypsin inhibitors is reduced. Cooking in boiling water and the use of microwaves have been shown to be effective. Sufficient pre-soaking of the seeds (~24 hours) is important, as the effect is greater in seeds with a higher moisture content.

Objective

Determine trypsin inhibitor activity in yellow and red soy.


Determine the effect of thermal treatment of soy (cooking in boiling water for 3 h (A); 3 h in ventilator oven at 104 °C (B); 6 h in ventilator oven at 104 °C (C)) on trypsin inhibitor activity.

Principle of method

The determination of trypsin inhibitor activity is based on the fact that trypsin catalyses the cleavage of *p*-nitroaniline from the synthetic substrate N-benzoyl-D,L-arginine-*p*-nitroanilide (BAPNA).

The trypsin inhibitor in the soybean extract reduces the activity of trypsin so that less product *p*-nitroaniline is formed.

How much less product is formed in a given time is proportional to the activity of the trypsin inhibitor or the content of the trypsin inhibitor in the reaction mixture.

N-benzoyl-D,L-arginine-*p*-nitroanilide (BAPNA)

N-benzoyl-D,L-arginine

p-nitroaniline

The content of the *p*-nitroaniline reaction product is determined spectrophotometrically.


Experimental procedure

- a) Extraction of trypsin inhibitor from soy (we perform the procedure of determining trypsin inhibitor activity in yellow and red soy, and in yellow soy that has been thermally treated in three different ways)
- b) Spectrophotometric determination of trypsin inhibitor activity

a) Extraction of trypsin inhibitor from soy (extract preparation)


- Weigh 1.0 g of ground and defatted soy ($m_{\text{defatted soy}}$) into a 100 mL beaker; soy is defatted with diethyl ether.
- Add 50 mL of the extraction solvent (NaOH solution (0.01 mol/L)); adjust the pH (HCl and NaOH solutions) to 8.6 (pH meter).
- Stirring on a magnetic stirrer for 30 min at room temperature; the trypsin inhibitor is extracted from the soy.
- Filter through a funnel with gauze into a measuring cylinder; the filtrate (trypsin inhibitor extract) is turbid; we need a clear extract, as the determination of activity is done spectrophotometrically; pour the filtrate (extract) into two centrifuges; make sure that the two centrifuges (together with their contents) are of equal weight; centrifugation follows.
- The supernatant (the clear part above the sediment after centrifugation in a centrifuge) is carefully poured into the measuring cylinder and the volume of the crude extract is read ($V_{\text{crude extr}}$).
- Prepare a diluted extract: measure a certain volume of crude extract into each individual 50 mL flask ($V_{\text{p crude extr}} = 1.0 \text{ mL}; 1.5 \text{ mL}; 2.0 \text{ mL}$), fill to the mark with distilled water, close the flasks and mix the contents; the volume of diluted extract ($V_{\text{diluted extr}}$) is 50 mL.

SCHEME OF TRYPSIN INHIBITOR EXTRACTION PROCEDURE

DILUTED EXTRACT PREPARATION SCHEME

pipette out crude extract into a 50 mL flask
($V_{p \text{ crude extr}} = 1.0 \text{ mL}; 1.5 \text{ mL}; 2.0 \text{ mL}$)

Different $V_{p \text{ crude extr}}$ mean different amounts of trypsin inhibitor in the reaction mixture in spectrophotometric determination.

b) Spectrophotometric determination of trypsin inhibitor activity


- Prepare three test tubes for each individual diluted extract: EP1 (blank sample), EP2 (without trypsin inh.), EP3 (contains trypsin inh.).
- Measure into a test tube EP1 (blank sample): 2.0 mL of the diluted extract, 1.0 mL of acetic acid solution (30%), 5.0 mL of BAPNA solution and 2.0 mL of trypsin solution (why this specified order?); mix with a vortex mixer.
- Measure a certain volume of diluted extract ($V_{p \text{ diluted extr}}$), water and trypsin (Table 1 on p. 154) into test tubes EP2 and EP3; mix with a vortex mixer.

Table 1: Volume of diluted extract ($V_{p \text{ diluted extr}}$), water ($V_{p \text{ water}}$) and trypsin ($V_{p \text{ trypsin}}$) for transfer to EP2 in EP3.

test tube	$V_{p \text{ diluted extr}}$ (mL)	$V_{p \text{ water}}$ (mL)	$V_{p \text{ trypsin}}$ (mL)
EP2	-	2.0	2.0
EP3	2.0	-	2.0

- The BAPNA solution and EP1, EP2 and EP3 with the described contents are placed in the water bath at 37 °C for approx. 5 min.
- Then pipette 5.0 mL of BAPNA solution (substrate) into EP2 and EP3.
- The reaction in EP2 and EP3 is allowed to run for exactly 10 min (water bath at 37 °C).
- Stop the reaction in EP2 and EP3 by adding 1.0 mL of acetic acid solution (30%) (why acetic acid solution?); mix with a vortex mixer.
- Pour the contents of the test tubes into cuvettes and measure the absorbance of the solution in EP2 (A_2) and EP3 (A_3) against a blank sample (EP1) at 410 nm.

SCHEME OF SPECTROPHOTOMETRIC DETERMINATION

Concentration of the product *p*-nitroaniline in EP2 > concentration of the product *p*-nitroaniline in EP3;
therefore: $A_2 > A_3$

Calculations

Trypsin inhibitor activity is proportional to ΔA ($\Delta A = A_2 - A_3$).

The greater the inhibition of trypsin, the less product is produced and the lower is A_3 , the greater is ΔA .

Trypsin inhibitor activity is expressed as the number of trypsin inhibition units (TIU).

By definition: 1 trypsin inhibition unit (1 TIU) corresponds to $\Delta A = 0.010$, determined under the conditions described.

$$\text{tryps. inh. activity in reaction mixture} = \frac{\Delta A}{0.010}$$

According to the experiment description (Table 1 on p. 154): a certain $V_{p \text{ diluted extr}}$ was measured into the reaction mixture (in EP3), therefore:

$$\text{tryps. inh. activity in } V_{p \text{ diluted extr}} = \text{tryps. inh. activity in reaction mixture}$$

$$\text{tryps. inh. activity / mL of diluted extr.} = \frac{\text{tryps. inh. activity in } V_{p \text{ diluted extr}}}{V_{p \text{ diluted extr}}}$$

tryps. inh. activity in total $V_{\text{diluted extr}} = V_{\text{diluted extr}} \times \text{tryps. inh. activity / mL of diluted extr.}$

According to the experiment description (p. 151, 152): $V_{\text{diluted extr}} = 50 \text{ mL}$; the diluted extract is prepared by measuring a certain $V_{\text{p crude extr}}$ and adding water, therefore:

tryps. inh. activity in $V_{\text{p crude extr}} = \text{tryps. inh. activity in total } V_{\text{diluted extr}}$

tryps. inh. activity / mL of crude extr. = $\frac{\text{tryps. inh. activity in } V_{\text{p crude extr}}}{V_{\text{p crude extr}}}$ (unit: TIU/mL)

According to the experiment description (p. 151, 153): the trypsin inhibitor was extracted from soybeans and after filtration and centrifugation, the $V_{\text{crude extr}}$ was read, therefore:

tryps. inh. activity in total $V_{\text{crude extr}} = V_{\text{crude extr}} \times \text{tryps. inh. activity / mL of crude extr.}$

tryps. inh. activity in defatted soy = tryps. inh. activity in total $V_{\text{crude extr}}$

tryps. inh. activity / mg of defatted soy = $\frac{\text{tryps. inh. activity in defatted soy}}{m_{\text{defatted soy}}}$ (unit: TIU/mg)

Example: an extract ($V_{\text{crude extr}} = 42 \text{ mL}$) was prepared from defatted soybeans ($m_{\text{defatted soy}} = 1000 \text{ mg}$) and according to the experiment procedure described ($V_{\text{p crude extr}} = 1.5 \text{ mL}$, $V_{\text{diluted extr}} = 50 \text{ mL}$, $V_{\text{p diluted extr}} = 2.0 \text{ mL}$) $\Delta A = 0.050$ was determined.

$$\text{tryps. inh. activity in reaction mixture} = \frac{\Delta A}{0.010} = \frac{0.050}{0.010} = 5.0 \text{ TIU}$$

$$\text{tryps. inh. activity in } V_{\text{p diluted extr}} = \text{tryps. inh. activity in reaction mixture} = 5.0 \text{ TIU}$$

$$\text{tryps. inh. activity / mL of diluted extr.} = \frac{\text{tryps. inh. activity in } V_{\text{p diluted extr}}}{V_{\text{p diluted extr}}} = \frac{5.0 \text{ TIU}}{2.0 \text{ mL}} = 2.5 \text{ TIU/mL}$$

$$\begin{aligned} \text{tryps. inh. activity in total } V_{\text{diluted extr}} &= V_{\text{diluted extr}} \times \text{tryps. inh. activity / mL of diluted extr.} = \\ &= 50 \text{ mL} \times 2.5 \text{ TIU/mL} = 125 \text{ TIU} \end{aligned}$$

$$\text{tryps. inh. activity in } V_{\text{p crude extr}} = \text{tryps. inh. activity in total } V_{\text{diluted extr}} = 125 \text{ TIU}$$

$$\begin{aligned} \text{tryps. inh. activity / mL of crude extr.} &= \frac{\text{tryps. inh. activity in } V_{p \text{ crude extr}}}{V_{p \text{ crude extr}}} = \frac{125 \text{ TIU}}{1.5 \text{ mL}} = \\ &= 83.3 \text{ TIU/mL} \end{aligned}$$

$$\begin{aligned} \text{tryps. inh. activity in total } V_{\text{crude extr}} &= V_{\text{crude extr}} \times \text{tryps. inh. activity / mL of crude extr.} = \\ &= 42 \text{ mL} \times 83.3 \text{ TIU/mL} = 3500 \text{ TIU} \end{aligned}$$

$$\text{tryps. inh. activity in defatted soy} = \text{tryps. inh. activity in total } V_{\text{crude extr}} = 3500 \text{ TIU}$$

$$\begin{aligned} \text{tryps. inh. activity / mg of defatted soy} &= \frac{\text{tryps. inh. activity in defatted soy}}{m_{\text{defatted soy}}} = \\ &= \frac{3500 \text{ TIU}}{1000 \text{ mg}} = 3.5 \text{ TIU/mg} \end{aligned}$$

Report

Experimental data

soy sample investigated:

$m_{\text{defatted soy}} = \dots\dots\dots$ $V_{\text{crude extr}} = \dots\dots\dots$ $V_{\text{p crude extr}} = \dots\dots\dots$ $A_2 = \dots\dots\dots$ $A_3 = \dots\dots\dots$

Further data can be found in the description of the experimental procedure.

Calculations

$\Delta A =$

tryps. inh. activity in reaction mixture =

tryps. inh. activity in $V_{\text{p diluted extr}} =$

tryps. inh. activity / mL of diluted extr. =

tryps. inh. activity in total $V_{\text{diluted extr}}$ =

tryps. inh. activity in $V_{\text{p crude extr}}$ =

tryps. inh. activity / mL of crude extr. =

tryps. inh. activity in total $V_{\text{crude extr}}$ =

tryps. inh. activity in defatted soy =

tryps. inh. activity / mg of defatted soy =

Complete the table and comment on the trypsin inhibitor activity in soy depending on the variety and thermal treatment.

Table 2: Trypsin inhibitor activity in soy.

	red soy			yellow soy			therm. treat. yellow soy (A)			therm. treat. yellow soy (B)			therm. treat. yellow soy (C)		
$V_{p \text{ crude extr}}$ (mL)	1.0	1.5	2.0	1.0	1.5	2.0	1.0	1.5	2.0	1.0	1.5	2.0	1.0	1.5	2.0
ΔA															
tryps. inh. activity in reaction mixture (TIU)															
tryps. inh. activity in $V_{p \text{ diluted extr}}$ (TIU)															
tryps. inh. activity / mL of diluted extr. (TIU/mL)															
tryps. inh. activity in total $V_{\text{diluted extr}}$ (TIU)															
tryps. inh. activity in $V_{p \text{ crude extr}}$ (TIU)															
tryps. inh. activity / mL of crude extr. (TIU/mL)															
tryps. inh. activity in total $V_{\text{crude extr}}$ (TIU)															
tryps. inh. activity in defatted soy (TIU)															
tryps. inh. activity / mg of defatted soy (TIU/mg)															

Comment

QUANTITATIVE DETERMINATION OF STARCH IN FLOUR

Starch

Starch is a mixture of two different polymers: amylose and amylopectin.

Amylose consists of several 100 to several 1000 glucose units linked by an α -(1→4)-glycosidic bond. The amylose molecule has one reducing end. The polymer forms a helix. Hydrogen bonds are formed between –OH groups of the glucose units in the polymer, which stabilize the helical structure. The molar mass of amylose depends on the botanical source (species, varieties). Amylose in wheat starch contains on average 1000-2000 glucose units and in potato starch up to 4500.

The amylopectin polymer is larger than amylose and consists of 60000 to 3000000 glucose units. The polymer is branched. A branch occurs approximately every 10 to 100 glucose units in the chain. In the polymer, the glucose units are linked by an α -(1→4)-glycosidic bond. At the branch point, the glucose molecule from C1 is linked to the main chain by an α -(1→6)-glycosidic bond. The side chains of amylopectin form double helices. The amylopectin molecule has one reducing end in the main chain.

Plants usually contain more amylopectin. The weight ratio between amylopectin and amylose depends on the botanical source and is generally 3:1.

Starch synthesis takes place in the organelles of the plant cells, in the chloroplasts of the leaves and in the amyloplasts, where it is stored in the form of starch granules. This so-called reserve starch, which is mainly found in the endosperm of seeds, and in roots and tubers, represents an energy store. The size and shape of the starch granules depend on the botanical source (species, varieties). The diameter of the starch granules in corn is 15 μm , in wheat 25 μm and in potatoes 40 μm .

Due to the close arrangement of the starch polymers in the granules, hydration in cold water is quite limited, which results in poor solubility. The structure of the granules does not change significantly. In cold water, the starch forms a suspension with low viscosity. Noticeable changes occur at higher temperatures (55-70 $^{\circ}\text{C}$), which weakens the strength of the hydrogen bonds between the polymer segments. This allows the polymers to separate from each other. The structure becomes more amorphous. As a result, more and more water penetrates the granules, which leads to increased hydration and swelling. Amylose molecules, and to a lesser extent, amylopectin molecules are released from the granules and the granules disintegrate.

The main sources of starch in the human diet are tubers, cereals, legumes, root vegetables, chestnuts, nuts and fruit. As the fruit ripens, the starch content decreases. Ripe fruit (with a few exceptions, such as bananas) contain no starch.

Starch is an important nutrient in the human diet. The digestion of starch begins in the mouth in the presence of the enzyme α -amylase in saliva, continues in the small intestine in the presence of α -amylase from the pancreas and then under the influence of maltase and isomaltase, which break down the resulting oligosaccharides and maltose into glucose.

Starch is an important industrial raw material. Modified starch is obtained from native starch which is chemically, physically and/or enzymatically modified to improve its functionality (more suitable technological and sensory properties). Native and modified starches are used in the food industry as thickeners, stabilizers of foams, emulsions and suspensions, as emulsifiers, gelling agents, humectants, fillers, cryostabilizers, encapsulation agents, carriers of aromas and fat substitutes.

Starch hydrolysates are products of starch modification, where glycosidic bonds between glucose units are cleaved by acid (HCl) and suitable enzymes. The products of hydrolysis are various dextrans, oligosaccharides, maltose and glucose. Acid hydrolysis is less specific than enzymatic hydrolysis. At industrial level, endoamylases (e.g. α -amylase) and exoamylases (e.g. β -amylase and amyloglucosidase) are used for the enzymatic hydrolysis of starch.

Objective

Determine the mass fraction of starch (w_{starch}) in wheat flour.


Determine the influence of the starch degradation procedure on the glucose content in the starch hydrolysate.

Principle of method

Soaking starch in boiling water enables better hydration and swelling of the granules, making them more susceptible to acids and enzymes and allowing more efficient starch degradation.

The end product of starch degradation by acid hydrolysis (it is less specific) and hydrolysis by the enzymes α -amylase (cleaves α -(1 \rightarrow 4)-glycosidic bonds within the chain) and amyloglucosidase (cleaves α -(1 \rightarrow 4)-glycosidic bonds from the non-reducing end of oligosaccharides and α -(1 \rightarrow 6)-glycosidic bonds (more slowly)) is glucose.

The determination of glucose in the hydrolysate is based on the fact that glucose is a strong reducing agent and reduces the nitro group in 3,5-dinitrosalicylic acid to the amino group.

The content of the reaction product 3-amino-5-nitrosalicylic acid is determined spectrophotometrically. The absorbance is proportional to the concentration of the reaction product. The concentration of the reaction product is proportional to the glucose content in the starch hydrolysate.

To evaluate how the starch degradation procedure (we omit certain steps) affects the glucose content in the hydrolysate, we perform the experiment in 8 different ways (Table 1 on p. 171).

Table 1: Starch degradation procedure.

	acid	heating	enzymes
1A	✓	✓	✓
1B	✓	✓	x
2A	✓	x	✓
2B	✓	x	x
3A	x	✓	✓
3B	x	✓	x
4A	x	x	✓
4B	x	x	x


Experimental procedure

- a) Acid hydrolysis
- b) Enzymatic hydrolysis
- c) Spectrophotometric determination of glucose in starch hydrolysate
- d) Calibration curve for determination of glucose in starch hydrolysate

a) Acid hydrolysis

- Weigh 0.80 g of flour (m_{flour}) in a 100 mL beaker.
- Add 25 mL of HCl solution (2 mol/L) and 25 mL of distilled water.
- Mixing with a glass rod at room temperature for 20 min.
- Place the beaker on a hot plate and wait for the solution to boil; boiling for 20 min.
- Cool to room temperature.
- Adjust the pH of the solution to 6.9 with NaOH solution (3 mol/L).
- Transfer the solution from the beaker quantitatively to a 100 mL volumetric flask, fill to the mark with distilled water, close the flask and mix the contents.


PROCEDURE SCHEME

b) Enzymatic hydrolysis

- Transfer 0.40 mL of the solution ($V_{p \text{ solution}}$) from the 100 mL flask to a test tube.
- Add 0.70 mL of the α -amylase solution.
- Leave for exactly 3 min.
- Add 0.10 mL of the amyloglucosidase solution.
- Incubation for 10 min in a water bath at 55 °C.
- Cool to room temperature.

SCHEME OF ENZYMATIC HYDROLYSIS PROCEDURE

Influence of starch degradation procedure on glucose content in hydrolysate:

1A and 1B: weight, add 25 mL acid, 25 mL water, room temp./20 min, boiling/20 min, cooling, neutralization, quantitatively to 100 mL volumetric flask, then:

1A: 0.40 mL into test tube, add 0.70 mL α -amylase, 0.10 mL amyloglucosidase, water bath

1B: 0.40 mL into test tube, add 0.80 mL water, water bath.

2A and 2B: weight, add 25 mL acid, 25 mL water, room temp./40 min, no cooling, neutralization, quantitatively to 100 mL volumetric flask, then:

2A: 0.40 mL into test tube, add 0.70 mL α -amylase, 0.10 mL amyloglucosidase, water bath

2B: 0.40 mL into test tube, add 0.80 mL water, water bath.

3A and 3B: weight, add 50 mL water, room temp./20 min, boiling/20 min, cooling, no neutralization, quantitatively to 100 mL volumetric flask, then:

3A: 0.40 mL into test tube, add 0.70 mL α -amylase, 0.10 mL amyloglucosidase, water bath

3B: 0.40 mL into test tube, add 0.80 mL water, water bath.

4A and 4B: weight, add 50 mL water, room temp./40 min, no cooling, no neutralization, quantitatively to 100 mL volumetric flask, then:


4A: 0.40 mL into test tube, add 0.70 mL α -amylase, 0.10 mL amyloglucosidase, water bath

3B: 0.40 mL into test tube, add 0.80 mL water, water bath.

c) Spectrophotometric determination of glucose in starch hydrolysate

- After complete hydrolysis, add 1.0 mL of 3,5-dinitrosalicylic acid solution (1 g/100 mL) to the starch hydrolysate in a test tube and mix.
- Prepare a blank sample in another test tube: pipette 1.2 mL of distilled water and 1.0 mL of 3,5-dinitrosalicylic acid solution into a test tube, mix.
- Place both test tubes in a boiling water bath for 5 min, then cool to room temperature.
- Transfer the contents of each test tube quantitatively to a 25 mL volumetric flask, fill to the mark with distilled water, close and mix the contents.
- Prepare two cuvettes: one cuvette for the sample (CU_{sample}) and one cuvette for the blank sample (CU_{blank}).
- Pour the contents of the flasks into a suitable cuvette.
- Measure the absorbance of the sample against a blank sample at 540 nm ($A_{\text{samp}540}$).

SCHEME OF SPECTROPHOTOMETRIC DETERMINATION OF GLUCOSE IN STARCH HYDROLYSATE

d) Calibration curve for determination of glucose in starch hydrolysate


- Use a glucose stock solution with a mass concentration of 2.0 mg/mL ($\gamma_{\text{glu stock}}$).
- Prepare 6 test tubes; pipette the appropriate volume of glucose stock solution ($V_{\text{p stock glu}}$) and distilled water ($V_{\text{p water}}$) into each test tube as indicated in Table 2 on p. 186.
- Add 1.0 mL of 3,5-dinitrosalicylic acid solution to each test tube, mix, place the test tubes in a boiling water bath for 5 min; then cool to room temperature.
- Prepare six 25 mL volumetric flasks; transfer the contents of each test tube quantitatively to the flask, fill to the mark with distilled water, close the flask and mix the contents.
- Prepare 6 cuvettes; pour the contents of each flask into a cuvette and measure A_{540} of all solutions against a blank sample; enter the data in Table 2 on p. 186 (m_{glucose} in 25 mL flask is calculated from $\gamma_{\text{glu stock}}$).

SCHEME OF CALIBRATION CURVE PREPARATION

prepare 6 test tubes
into each:

$V_{p \text{ stock glu}}$ ($V_{\text{glu stock}} = 2.0 \text{ mg/mL}$)
 $V_{p \text{ water}}$
1.0 mL 3,5-dinitrosal. acid


blank sample:
1.2 mL water
1.0 mL 3,5-dinitrosal. acid

boiling water/5 min

cool

quantitatively
transfer to 25 mL
volumetric flask

Calculations


Taking into account the values in Table 2 on p. 186, we draw a graph of the dependence of A_{540} on m_{glucose} in 25 mL flask. Read from the graph the corresponding value of m_{glucose} in 25 mL flask at measured $A_{\text{samp}540}$ for the investigated starch hydrolysate.

For example, if we determine the value $A_{\text{samp}540} = 0.316$ for the starch hydrolysate, we read m_{glucose} in 25 mL flask as shown on p. 183.

The dependence of A_{540} on m_{glucose} in 25 mL flask can also be described by the equation of a straight line: $A_{540} = k \times m_{\text{glucose}}$ in 25 mL flask.

The slope (k) is determined by linear regression analysis of the values of A_{540} and m_{glucose} in 25 mL flask (Table 2 on p. 186). For the investigated starch hydrolysate, we calculate the corresponding value of m_{glucose} in 25 mL flask at measured $A_{\text{samp}540}$. The linear regression analysis is carried out using appropriate computer programs (MS Excel or Origin).

How to calculate the mass of starch in the weight of flour (m_{starch}) from the m_{glucose} in the starch hydrolysate in 25 mL flask is explained on p. 184, 185.

Dependence of A_{540} on m_{glucose} in 25 mL flask

m_{glucose} in starch hydrolysate in 25 mL flask = m_{glucose} in starch hydrolysate in test tube

m_{glucose} in starch hydrolysate in test tube = m_{glucose} in $V_{\text{p solution}}$

According to the experiment description: 0.40 mL of solution is measured from 100 mL flask ($V_{\text{p solution}} = 0.40$ mL) into a test tube, therefore:

$$m_{\text{glucose}} \text{ in 100 mL flask} = \frac{m_{\text{glucose}} \text{ in } V_{\text{p solution}}}{V_{\text{p solution}}} \times 100 \text{ mL}$$

We take into account the fact that water is released when a glycosidic bond is formed between glucose units in the polysaccharide ($M_{\text{glucose}} = 180$ g/mol, $M_{\text{water}} = 18$ g/mol), therefore:

$M_{\text{glucose unit}}$ in polysaccharide = $M_{\text{glucose}} \times 0.9$ or: $m_{\text{glucose units}}$ in starch polymer = $m_{\text{glucose}} \times 0.9$

$m_{\text{glucose units}}$ in 100 mL flask = m_{glucose} in 100 mL flask $\times 0.9$

According to the experiment description: when preparing the hydrolysate, the contents of the beaker are transferred quantitatively to a 100 mL volumetric flask, therefore:

$$m_{\text{glucose units}} \text{ in 100 mL flask} = m_{\text{glucose units}} \text{ in starch polymer in flour} = \\ = m_{\text{starch}} \text{ in flour}$$

$$w_{\text{starch}} \text{ in flour} = \frac{m_{\text{starch}} \text{ in flour}}{m_{\text{flour}}} \times 100 \%$$

Report

Analysis of starch hydrolysate

starch degradation process:


$m_{\text{flour}} = \dots\dots\dots A_{\text{samp540}} = \dots\dots\dots$

Further data can be found in the description of the experimental procedure.

Calibration curve

Table 2: Dependence of A_{540} on m_{glucose} in 25 mL flask.

test tube	$V_{\text{p stock glu}}$ (mL)	m_{glucose} in 25 mL flask (mg)	$V_{\text{p water}}$ (mL)	A_{540}
1	0.20	0.40	1.00	
2	0.40	0.80	0.80	
3	0.60	1.20	0.60	
4	0.80	1.60	0.40	
5	1.00	2.00	0.20	
6	1.20	2.40	0.00	

Dependence of A_{540} on m_{glucose} in 25 mL flask

Calculation of m_{glucose} in starch hydrolysate and the mass fraction of starch in the flour (w_{starch})

The value of m_{glucose} in the starch hydrolysate in 25 mL flask read from the calibration curve or calculated using k ($k = \dots\dots\dots$) at measured A_{samp540} : $\dots\dots\dots$

Calculate w_{starch} ; write the complete calculation with numbers and units.

m_{glucose} in starch hydrolysate in test tube =

m_{glucose} in 100 mL flask =

$m_{\text{glucose units}}$ in starch polymer =

m_{starch} in flour =

w_{starch} in flour =

Complete the table and comment on the effect of the starch degradation procedure on the glucose content in the starch hydrolysate and the calculated mass fraction of starch.

Table 3: Glucose content in prepared starch hydrolysates and calculated mass fraction of starch (w_{starch}).

starch degradation procedure	A_{samp540}	m_{glucose} in starch hydrolysate in 25 mL flask (mg)	w_{starch} (%)
1A			
1B			
2A			
2B			
3A			
3B			
4A			
4B			

Comment

QUANTITATIVE DETERMINATION OF CHOLESTEROL IN MILK

Cholesterol

Cholesterol belongs to a group of lipids called steroids. Sterols are compounds with an –OH group attached to the steroid skeleton.

Cholesterol is a component of cell membranes and influences the fluidity or stiffness of the membranes. It is present in large quantities in nervous tissue, especially in the brain. It is also a precursor for the synthesis of bile acids, vitamin D3 (cholecalciferol) and steroid hormones. Humans can synthesize cholesterol themselves (mainly in the liver) or consume it with food. Due to its poor water solubility, it is transported through the body in the blood and other extracellular fluids as part of lipid complexes with proteins, i.e. lipoproteins.

Cholesterol belongs to the so-called animal sterols and is found in foods such as meat, fat, egg yolk, milk and dairy products. In milk, cholesterol is found in lipoprotein complexes, mainly in the phospholipid monolayer that surrounds the fat globules, and to a lesser extent (as esters with fatty acids) in their interior. The cholesterol content in milk is related to the milk fat content. If dietary intake of cholesterol is too high or metabolic disorders occur, the concentration of cholesterol in the blood rises, which represents an increased risk of cardiovascular disease. This underlines the importance of determining cholesterol in food.

Plants synthesize various so-called plant sterols, phytosterols (β -sitosterol, stigmasterol, ...). They are mainly found in vegetable oils. When ingested, phytosterols are supposed to reduce the absorption of cholesterol in the digestive tract and thus lower the concentration of cholesterol in the blood. For this reason, they are added to some foods.

Objective

Determine the mass concentration of cholesterol in milk samples (Y_{cholM}) with 0.5%, 1.5% or 3.5% milk fat.

Find out how the cholesterol concentration in milk is related to the content of milk fat.

Principle of method

In the milk sample, the triacylglycerols are degraded by saponification and the cholesterol is released from cholesterol esters and the lipoprotein complex.

After the free cholesterol has been extracted into the organic phase, the product formed in the following reaction in the presence of a strong acid is determined spectrophotometrically:

The absorbance is proportional to the concentration of the reaction product and the concentration of the reaction product is proportional to the cholesterol content.

Experimental procedure

- a) Saponification with KOH in ethanol
- b) Extraction of cholesterol into hexane and hexane evaporation
- c) Spectrophotometric determination of cholesterol
- d) Calibration curve for cholesterol determination

a) Saponification with KOH in ethanol

- Measure 2 mL of milk ($V_{p, \text{milk}}$) into a glass-stoppered test tube, add 3 mL of ethanol (95%) and 2 mL of KOH solution (50%); close the test tube and mix the contents.
- Heat for 15 min at 60 °C in a water bath; the triacylglycerols are degraded by saponification and the cholesterol is released from the lipoprotein complex.
- Cool to room temperature.


b) Extraction of cholesterol into hexane and hexane evaporation

- Add 5 mL of hexane ($V_{\text{hex ph}}$) and 3 mL of water.
- Close the test tube and shake it; the cholesterol is extracted from the aqueous phase into the hexane.
- Then place the test tube in a test tube rack and wait until the aqueous and hexane phases are separated; the aqueous phase is in the lower part of the test tube, the hexane phase is above it; the interface between the two phases is visible.
- Transfer the hexane phase to the next clean test tube using a dropper; make sure that you do not catch the aqueous phase with the dropper.
- Transfer the appropriate volume of the hexane phase ($V_{\text{p hex ph}}$) to the next clean test tube using a measuring pipette; for milk with 0.5% fat, $V_{\text{p hex ph}} = 3 \text{ mL}$; for milk with 1.5% or 3.5% fat, $V_{\text{p hex ph}} = 2 \text{ mL}$.
- Remove the hexane by evaporation in a vacuum dryer (30 min, 40 °C, 30 mbar).

c) Spectrophotometric determination of cholesterol

- To the dry residue after hexane evaporation (contains cholesterol), add 4 mL of *o*-phthalaldehyde solution in concentrated acetic acid, mix carefully; after 10 min, add 2 mL of concentrated H₂SO₄ solution (be careful when working!).
- After 30 min, measure the absorbance at 550 nm ($A_{\text{samp}550}$) against a blank sample.
- The blank sample consists of 4 mL of *o*-phthalaldehyde solution and 2 mL of concentrated H₂SO₄ solution.


PROCEDURE SCHEME

d) Calibration curve for cholesterol determination

- Use a cholesterol stock solution in 95% ethanol with a mass concentration of 0.10 mg/mL ($Y_{\text{chol stock}}$).
- Prepare 5 glass-stoppered test tubes; pipette the appropriate volume of cholesterol stock solution ($V_{\text{p chol stock}}$) and 95% ethanol ($V_{\text{95% ethanol}}$) into each test tube as indicated in Table 1 on p. 206.
- Add 2 mL of KOH solution (50%) to each test tube, heat for 15 min at 60 °C in a water bath, and proceed similarly to the determination of cholesterol in milk.
- The only difference is that we transfer 4 mL of the hexane phase ($V_{\text{p hex ph}} = 4 \text{ mL}$) before the hexane evaporation and not 2 mL or 3 mL as in the determination of cholesterol in milk.
- Measure A_{550} of all solutions and enter the data in Table 1 on p. 206; how to calculate the cholesterol concentration in the reaction mixture ($Y_{\text{chol react mixt}}$) for the calibration curve is explained on p. 203.

SCHEME OF CALIBRATION CURVE PREPARATION

Calculations


Taking into account the values in Table 1 on p. 206, we draw a graph of the dependence of A_{550} on $Y_{\text{chol react mixt}}$. Read from the graph the corresponding value of $Y_{\text{chol react mixt}}$ for the milk sample at measured A_{samp550} .

For example, if we determine the value $A_{\text{samp550}} = 0.525$, we read the value of $Y_{\text{chol react mixt}}$ as shown on p. 202.

The dependence of A_{550} on $Y_{\text{chol react mixt}}$ can also be described by the equation of a straight line:
$$A_{550} = k \times Y_{\text{chol react mixt}}$$

The slope (k) is determined by linear regression analysis of the values of A_{550} and $Y_{\text{chol react mixt}}$ (Table 1 on p. 206). We calculate the corresponding value of $Y_{\text{chol react mixt}}$ for the milk sample at measured A_{samp550} . The linear regression analysis is carried out using appropriate computer programs (MS Excel or Origin).

How to calculate Y_{cholIM} for the investigated milk sample from the $Y_{\text{chol react mixt}}$ is explained on p. 204, 205.

Dependence of A_{550} on $Y_{\text{chol react mixt}}$

Calculation of cholesterol concentration in the reaction mixture ($Y_{\text{chol react mixt}}$) for the calibration curve

$$m_{\text{chol}} \text{ in } V_{\text{p chol stock}} = Y_{\text{chol stock}} \times V_{\text{p chol stock}}$$

Values for $V_{\text{p chol stock}}$ are presented in Table 1 on p. 206; $Y_{\text{chol stock}} = 0.10 \text{ mg/mL}$.

Cholesterol is extracted from the stock solution into hexane, therefore:

$$m_{\text{chol}} \text{ in total volume of hexane phase} = m_{\text{chol}} \text{ in } V_{\text{p chol stock}}$$

$$Y_{\text{chol hex ph}} = \frac{m_{\text{chol}} \text{ in total volume of hexane phase}}{V_{\text{hex ph}}}$$

According to the experiment description: $V_{\text{hex ph}} = 5 \text{ mL}$; when preparing the calibration curve, $V_{\text{p hex ph}} = 4 \text{ mL}$ is transferred to a new test tube.

$$m_{\text{chol}} \text{ in } V_{\text{p hex ph}} = Y_{\text{chol hex ph}} \times V_{\text{p hex ph}}$$

$m_{\text{chol}} \text{ in } V_{\text{p hex ph}} = m_{\text{chol}} \text{ before hexane evaporation} = m_{\text{chol}} \text{ after hexane evaporation} = m_{\text{chol}} \text{ in reaction mixture}$, therefore:

$$Y_{\text{chol react mixt}} = \frac{m_{\text{chol}} \text{ in reaction mixture}}{V_{\text{react mixt}}}$$

According to the experiment description: $V_{\text{react mixt}} = 6 \text{ mL}$.

Calculation of cholesterol concentration in milk (y_{cholM})

Read from the calibration curve or use k to calculate the corresponding value of $y_{\text{chol react mixt}}$ for the investigated milk sample at the measured A_{samp550} .

$$m_{\text{chol}} \text{ in reaction mixture} = y_{\text{chol react mixt}} \times V_{\text{react mixt}}$$

According to the experiment description: $V_{\text{react mixt}} = 6 \text{ mL}$.

m_{chol} in reaction mixture = m_{chol} after hexane evaporation = m_{chol} before hexane evaporation = m_{chol} in $V_{\text{p hex ph}}$, therefore:

$$y_{\text{chol hex ph}} = \frac{m_{\text{chol}} \text{ in } V_{\text{p hex ph}}}{V_{\text{p hex ph}}}$$

According to the experiment description: for milk with 0.5% fat, $V_{\text{p hex ph}} = 3 \text{ mL}$; for milk with 1.5% or 3.5% fat, $V_{\text{p hex ph}} = 2 \text{ mL}$.

$$m_{\text{chol}} \text{ in total volume of hexane phase} = y_{\text{chol hex ph}} \times V_{\text{hex ph}}$$

According to the experiment description: $V_{\text{hex ph}} = 5 \text{ mL}$.

The cholesterol is extracted from the milk into hexane, therefore:

$$m_{\text{chol}} \text{ in } V_{\text{p milk}} = m_{\text{chol}} \text{ in total volume of hexane phase}$$

$$Y_{\text{cholM}} = \frac{m_{\text{chol}} \text{ in } V_{\text{p milk}}}{V_{\text{p milk}}}$$

Report

Milk sample analysis

milk fat content =

$V_{p \text{ hex ph}} = \dots\dots\dots A_{\text{samp550}} = \dots\dots\dots$

Further data can be found in the description of the experimental procedure.

Calibration curve

Table 1: Dependence of A_{550} on $Y_{\text{chol react mixt}}$ *

test tube	$V_{p \text{ chol stock}}$ (mL)	$V_{95\% \text{ ethanol}}$ (mL)	$Y_{\text{chol react mixt}}$ (mg/mL)	A_{550}
1	0.25	2.75	0.0033	
2	0.50	2.50	0.0067	
3	0.75	2.25	0.0100	
4	1.00	2.00	0.0133	
5	1.25	1.75	0.0167	

Calculation of cholesterol concentration in the reaction mixture ($y_{\text{chol react mixt}}$) for the calibration curve

Calculate the $y_{\text{chol react mixt}}$ for each of the tubes, and check whether your calculated value agrees with the value in Table 1 on p. 206.

$$m_{\text{chol}} \text{ in } V_{\text{p chol stock}} =$$


$$m_{\text{chol}} \text{ in total volume of hexane phase} =$$

$$Y_{\text{chol hex ph}} =$$

$$m_{\text{chol}} \text{ in } V_{\text{p hex ph}} =$$

$$m_{\text{chol}} \text{ in reaction mixture} =$$

$$Y_{\text{chol react mixt}} =$$

Dependence of A_{550} on $Y_{\text{chol react mixt}}$

Calculation of cholesterol concentration in milk (y_{cholM})

The value of $y_{\text{chol react mixt}}$ for the investigated milk sample read from the calibration curve or calculated using k ($k = \dots\dots\dots$) at measured A_{samp550} : $\dots\dots\dots$

Calculate y_{cholM} ; write the complete calculation with numbers and units.

m_{chol} in reaction mixture =

m_{chol} in $V_{\text{p hex ph}}$ =

$Y_{\text{chol hex ph}}$ =

m_{chol} in total volume of hexane phase =

m_{chol} in $V_{\text{p milk}}$ =

Y_{cholM} =

Complete the table and comment on the cholesterol concentration determined in relation to the milk fat content.

Table 2: Dependence of cholesterol concentration on milk fat content.

milk fat content (%)	0.5	1.5	3.5
$Y_{\text{chol react mixt}}$ (mg/mL)			
m_{chol} in reaction mixture (mg)			
m_{chol} in $V_{\text{p hex ph}}$ (mg)			
$Y_{\text{chol hex ph}}$ (mg/mL)			
m_{chol} in total volume of hexane phase (mg)			
m_{chol} in $V_{\text{p milk}}$ (mg)			
Y_{cholM} (mg/mL)			


Comment

DETERMINATION OF POLYPHENOL OXIDASE ACTIVITY IN POTATOES


Phenol oxidases

Phenol oxidases are enzymes that catalyse the aerobic oxidation of phenolic compounds, which leads to the formation of various brown-coloured high molecular weight products.

The first reaction is the hydroxylation of monophenols to *o*-diphenols:


the second reaction is the oxidation of *o*-diphenols to *o*-quinones:

Monophenols are compounds with one –OH group attached to a benzene ring. *o*-diphenols (catechol) are compounds where two –OH groups are attached to a benzene ring in the *ortho* position.

Enzymes that catalyse the hydroxylation of monophenols to *o*-diphenols belong to the monophenol oxidases. Enzymes that catalyse the oxidation of *o*-diphenols to *o*-quinones belong to the diphenol (catechol) oxidases. Enzymes that catalyse both, the hydroxylation of monophenols to *o*-diphenols and their oxidation to *o*-quinones belong to the tyrosinases.

Copper ions are involved in the catalytic activity of these enzymes.

The resulting *o*-quinones are reactive molecules. Brown-coloured products with high molecular weight are formed by the polymerization of *o*-quinones. Amines and proteins are involved in the polymerization reaction, which is not enzyme-catalysed.

Phenol oxidases are found in microorganisms as well as in plant and animal cells, where their activity leads to the formation of melanin pigments. The role of mentioned enzymes in plants is believed to be for protection against insects and pathogens. In a plant cell, they are found in the chloroplasts and chromoplasts. In the event of tissue damage, they are released from these organelles, come into contact with the substrate (phenolic compounds) and are activated by various factors. The resulting *o*-quinones can deactivate enzymes in an invading organism. Polymer products of *o*-quinones may also provide physical protection.

In the food industry, we pay attention to phenol oxidases because they cause enzymatic browning of fruit and vegetables during processing and storage. The browning occurs on the surface of the tissue that is exposed to mechanical processing or damage (peeling, cutting, grinding, ...). Browning can also occur in fruit or vegetable juice, puree, ... Due to the deterioration of sensory properties (appearance), such foods are less acceptable to the consumer.

The resulting waste causes economic losses and represents a burden on the environment. In some cases, such as raisins, prunes, cocoa, black tea and coffee, enzymatic browning is desirable.

Optimal conditions for enzyme activity are a pH range between 4.0 and 7.0 and a temperature between 30 and 50 °C. The enzyme is also relatively stable in the temperature range from 55 to 80 °C. It is therefore possible that the catalytic activity of phenol oxidases can also be demonstrated during a shorter thermal treatment, as the cell membranes are damaged at ~60 °C, so that the enzyme can come into contact with the substrate.

Enzymatic browning is mitigated/prevented by various physical and chemical treatments.

Physical treatments: dehydration, freezing, thermal treatment (blanching) to deactivate the enzyme; packaging in modified atmosphere, coating with sugar syrup, use of protective films to limit O₂ availability.

Chemical treatments to inhibit the enzymes include the use of acids (citric acid, malic acid, phosphoric(V) acid) to lower the pH, as phenol oxidases are weakly active at pH < 3; the addition of chelators (EDTA, oxalic acid, citric acid) that bind copper ions; the addition of reducing agents (sulfite, ascorbic acid, cysteine) that reduce *o*-quinones to diphenols.

Phenolic compounds

This is a large group of structurally diverse compounds. Phenolic compounds include: simple phenols or benzoquinones, phenolic acids, naphthoquinones, xanthenes, stilbenes, flavonoids, lignans, biflavonoids, lignins, coumarins and condensed tannins. Nutritionally, the most important are phenolic acids, flavonoids and tannins. Flavonoids are the largest group of phenolic compounds.

The distribution of these compounds in the cell/plant is not uniform. They are located in the vacuole (free, conjugated) and are bound to the cell wall. They accumulate in the outer layers and in the epidermal tissue of the plant. These compounds offer the plant protection against external stress factors, serve as visual markers and are important for growth and reproduction.

Phenolic compounds influence the sensory properties (colour, taste, aroma) of foods. Spices, berries, tea and olive oil are particularly rich in phenolic compounds.

Objective


Determine the polyphenol oxidase activity in old and new potatoes.

Determine the effect of thermal treatment of potatoes (cooking in boiling water for 45 min) on polyphenol oxidase activity.

Principle of method

The polyphenol oxidase (PPO) in potato extract catalyses the aerobic oxidation of the phenolic compound 3,4-dihydroxyphenylalanine (DOPA).

The rate of formation of the reaction product (*o*-quinone) is proportional to the activity of the polyphenol oxidase and to the content of polyphenol oxidase in the reaction mixture.

The concentration of the reaction product (o-quinone) and the rate of its formation are monitored spectrophotometrically.


Experimental procedure

- a) Extraction of polyphenol oxidase from potatoes (we perform the procedure of determining polyphenol oxidase activity in old and new potatoes and in thermally-treated new potatoes)
- b) Spectrophotometric determination of polyphenol oxidase activity

a) Extraction of polyphenol oxidase from potato (extract preparation)

- Peel the potatoes, cut them and weigh 25 g (m_{potato}) in a chopper.
- Add 50 mL of ice-cold phosphate buffer I (0.1 mol/L, pH 6.8; contains NaF).
- Homogenization for 1 min.
- Filter through a funnel with gauze into a measuring cylinder; the filtrate (PPO extract) is turbid; we need a clear extract, as the activity is determined spectrophotometrically; pour the filtrate (extract) into two centrifuges; make sure that the two centrifuges (together with their contents) are of equal weight; centrifugation follows.
- The supernatant (the clear part above the sediment after centrifugation in a centrifuge) is carefully poured into the measuring cylinder and the volume of the crude extract is read ($V_{\text{crude extr}}$).
- Prepare a diluted extract: measure 5 mL of the crude extract (that is $V_{\text{p crude extr}}$) into a test tube, add 5 mL of phosphate buffer II (0.1 mol/L, pH 6.8) and mix with a vortex mixer; the volume of the diluted extract ($V_{\text{diluted extr}}$) is 10 mL.

SCHEME OF POLYPHENOL OXIDASE EXTRACTION PROCEDURE

b) Spectrophotometric determination of polyphenol oxidase activity

- For each diluted extract prepare: sample cuvettes (CU_{sample}), control sample cuvette (CU_{contr}) and blank sample cuvette (CU_{blank}).
- The cuvette CU_{blank} contains water.
- Add buffer II and DOPA solution (Table 1 on p. 224) into the cuvette CU_{contr} , mix and measure the absorbance at 475 nm (A_{contr}) against a blank sample.
- Add buffer II, DOPA solution and a certain volume of diluted extract ($V_{\text{p diluted extr}}$) into cuvettes CU_{sample} (Table 1 on p. 224); the extract is added at the spectrophotometer shortly before starting the measurements.
- After adding the extract, start the stopwatch, mix the contents well, insert the cuvette into the spectrophotometer and start measuring the absorbance at 475 nm against a blank sample (A_{measured}).
- A_{measured} is monitored every 15 s until $t = 90$ s, then every 30 s until $t = 300$ s.

Table 1: Volume of buffer II ($V_{p \text{ buffer II}}$), DOPA solution ($V_{p \text{ DOPA}}$) and diluted extract ($V_{p \text{ diluted extr}}$) for transfer to cuvettes.

cuvette	$V_{p \text{ buffer II}}$ (mL)	$V_{p \text{ DOPA}}$ (mL)	$V_{p \text{ diluted extr}}$ (mL)
CU _{contr}	2.00	1.00	-
CU _{sample 1}	1.95	1.00	0.05
CU _{sample 2}	1.90	1.00	0.10
CU _{sample 3}	1.80	1.00	0.20
CU _{sample 4}	1.70	1.00	0.30

Different $V_{p \text{ diluted extr}}$ mean different amounts of enzyme in the reaction mixture.

Calculations

The absorbance of the reaction product (A_{sample}) is calculated by subtracting the A_{contr} value from the A_{measured} value ($A_{\text{sample}} = A_{\text{measured}} - A_{\text{contr}}$). Enter the data in Table 2 on p. 231 and draw a graph of the A_{sample} dependence on the incubation time.

To determine the initial rate, draw a tangent to the curve at time $t = 0$, as shown on p. 226, and read A_{sample} at 60 sec (1 min); this is $\Delta A_{\text{sample}} / \text{min}$.


PPO activity is proportional to the amount of a product in a given volume of the reaction mixture ($V_{\text{react mixt}}$) formed in one minute ($n_{\text{product}} / \text{min}$), and $n_{\text{product}} / \text{min}$ is proportional to $\Delta A_{\text{sample}} / \text{min}$.

$n_{\text{product}} / \text{min}$ is calculated taking into account the Beer-Lambert law ($A = \varepsilon \times l \times c$).

The concentration of the product in the reaction mixture (c_{product}) is given as:

$$c_{\text{product}} = \frac{n_{\text{product}}}{V_{\text{react mixt}}} = \frac{A_{\text{sample}}}{\varepsilon \times l} \quad \text{or:} \quad \frac{n_{\text{product}} / \text{min}}{V_{\text{react mixt}}} = \frac{\Delta A_{\text{sample}} / \text{min}}{\varepsilon \times l}$$

ε is molar absorption coefficient ($5012 \text{ L} \times \text{mol}^{-1} \times \text{cm}^{-1}$); l is the length of the optical path (dimension of the cuvette 1 cm).

Dependence of A_{sample} on incubation time

$$n_{\text{product}} / \text{min} = \frac{\Delta A_{\text{sample}} / \text{min}}{\varepsilon \times l} \times V_{\text{react mixt}}$$

The volume of the reaction mixture ($V_{\text{react mixt}}$) is the sum of the volumes pipetted into a cuvette $\text{CU}_{\text{sample}}$; according to the experiment description: $V_{\text{react mixt}} = 3 \text{ mL}$.

n_{product} express in μmol ($1 \mu\text{mol} = 1 \times 10^{-6} \text{ mol}$)

The polyphenol oxidase activity, which is proportional to $n_{\text{product}} / \text{min}$, is expressed as the number of enzyme units (U).

By definition: 1 unit of enzyme activity (1 U) catalyses the formation of 1 μmol of a product in 1 min, determined under the conditions described, therefore:

$$\text{PPO activity in reaction mixture} = \frac{1 \text{ U} \times n_{\text{product}} / \text{min}}{1 \mu\text{mol}/\text{min}}$$

For example, if we calculate that 0.02 μmol of product was formed in 1 min ($n_{\text{product}} / \text{min} = 0.02 \mu\text{mol}/\text{min}$), the PPO activity in the reaction mixture = 0.02 U.

According to the experiment description (Table 1 on p. 224): a certain $V_{p \text{ diluted extr}}$ was measured into a cuvette CU_{sample} (reaction mixture), therefore:

PPO activity in $V_{p \text{ diluted extr}}$ = PPO activity in reaction mixture

$$\text{PPO activity / mL of diluted extr.} = \frac{\text{PPO activity in } V_{p \text{ diluted extr}}}{V_{p \text{ diluted extr}}}$$

PPO activity in total $V_{\text{diluted extr}}$ = $V_{\text{diluted extr}} \times$ PPO activity / mL of diluted extr.

According to the experiment description (p. 221, 222): $V_{\text{diluted extr}} = 10 \text{ mL}$; diluted extract is prepared by measuring 5 mL of crude extract ($V_{p \text{ crude extr}}$) + 5 mL of buffer II, therefore:

PPO activity in $V_{p \text{ crude extr}}$ = PPO activity in total $V_{\text{diluted extr}}$

$$\text{PPO activity / mL of crude extr.} = \frac{\text{PPO activity in } V_{p \text{ crude extr}}}{V_{p \text{ crude extr}}} \quad (\text{unit: U/mL})$$

According to the experiment description (p. 221, 222): PPO was extracted from the potato and after filtration and centrifugation, the $V_{\text{crude extr}}$ was read, therefore:

$$\text{PPO activity in total } V_{\text{crude extr}} = V_{\text{crude extr}} \times \text{PPO activity / mL crude extr.}$$

$$\text{PPO activity in potato} = \text{PPO activity in total } V_{\text{crude extr}}$$

$$\text{PPO activity / mg of potato} = \frac{\text{PPO activity in potato}}{m_{\text{potato}}} \quad (\text{unit: U/mg})$$

Report

Experimental data

potato sample investigated:

$m_{\text{potato}} = \dots\dots\dots$

$V_{\text{crude extr}} = \dots\dots\dots$

$V_{\text{p crude extr}} = \dots\dots\dots$


$A_{\text{contr}} = \dots\dots\dots$

Further data can be found in the description of the experimental procedure.

Table 2: Dependence of A_{sample} on incubation time.

V_p diluted extr (mL)	old potato				new potato				therm. treated new potato			
	0.05	0.10	0.20	0.30	0.05	0.10	0.20	0.30	0.05	0.10	0.20	0.30
t(s)	A_{sample}	A_{sample}	A_{sample}	A_{sample}	A_{sample}	A_{sample}	A_{sample}	A_{sample}	A_{sample}	A_{sample}	A_{sample}	A_{sample}
0												
15												
30												
45												
60												
75												
90												
120												
150												
180												
210												
240												
270												
300												

Calculations and graph

Dependence of A_{sample} on incubation time

$\Delta A_{\text{sample}} / \text{min} =$

$n_{\text{product}} / \text{min} =$

PPO activity in reaction mixture =

PPO activity in $V_{\text{p diluted extr}}$ =

PPO activity / mL of diluted extr. =

PPO activity in total $V_{\text{diluted extr}}$ =

PPO activity in $V_{\text{p crude extr}}$ =

PPO activity / mL of crude extr. =

PPO activity in total $V_{\text{crude extr}}$ =

PPO activity in potato =

PPO activity / mg of potato =

Complete the table and comment on the PPO activity depending on the selected potato and thermal treatment.

Table 3: PPO activity in potato.

	old potato				new potato				therm. treated new potato			
$V_{p \text{ diluted extr}}$ (mL)	0.05	0.10	0.20	0.30	0.05	0.10	0.20	0.30	0.05	0.10	0.20	0.30
$\Delta A_{\text{sample}} / \text{min}$												
$n_{\text{product}} / \text{min}$ ($\mu\text{mol}/\text{min}$)												
PPO activity in reaction mixture (U)												
PPO activity in $V_{p \text{ diluted extr}}$ (U)												
PPO activity / mL of diluted extr. (U/mL)												
PPO activity in total $V_{\text{diluted extr}}$ (U)												
PPO activity in $V_{p \text{ crude extr}}$ (U)												
PPO activity / mL of crude extr. (U/mL)												
PPO activity in total $V_{\text{crude extr}}$ (U)												
PPO activity in potato (U)												
PPO activity / mg of potato (U/mg)												

Comment

REFERENCES

- Abram V., Zelenik-Blatnik M. 2002. Vaje iz živilske kemije za študente živilske tehnologije. 2nd ed. Ljubljana, Biotehniška fakulteta, Oddelek za živilstvo, <https://repozitorij.uni-lj.si/IzpisGradiva.php?id=122506>
- Fennema's Food Chemistry. 2008. Damodaran S., Parkin K. L., Fennema O. R. (eds.). 4th ed. Boca Raton, CRC Press, Taylor & Francis Group
- Pomeranz Y., Meloan C.E. 1980. Food Analysis. Westport, Connecticut, Avi Publishing Company
- Matissek R., Schnepel F.M., Steiner G. 1992. Lebensmittelanalytik. Springer-Verlag, Berlin/Heidelberg
- Miller D.D. 1998. Food Chemistry: A Laboratory Manual. New York [etc.], John Wiley & Sons

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© Univerza v Ljubljani, 2023. Vse pravice pridržane / University of Ljubljana, 2023. All rights reserved